

CTS-11030

Kit, PM, Pump



Technical Data Sheet

1.0 Description

This PM Pump kit contains all the parts necessary to keep Waters® 2690 and 2695 HPLC systems working at peak performance.

2.0 Specifications

2.1 Components

Qty	Part Description
2	Plunger, Sapphire
1	Seal, Plunger Wash, 2pk
1	Seal, Head Plunger, 2pk
1	Seal, Wash Tube, 4pk
1	Cartridge, Check Valve, 2pk
1	Filter, Insert
1	Battery
1	Sparge Diffuser
1	Seal, Face, 4pk
1	Kit, Seal Pack with Needle
1	Syringe, 250ul
4	Filter, Solvent, 10 Micron

2.2 Wetted Materials

UHMW-PE
Teflon®
316 Stainless Steel
Sapphire
Ruby
Borosilicate glass
Kel-F



3.0 Installation Information

CTS's parts are designed to be compatible and interchangeable with the OEM parts.

Please refer to your instrument service or maintenance manual in order to follow the exact installation procedure for this part.

4.0 Warranty

We view your business with us as a partnership based on mutual trust.

We extend a Lifetime Warranty against any manufacturing defects on all CTS parts. On the rare occasion that you experience a problem, we invite you to contact your Account Manager directly.

Returns can only be accepted if a Return Material Authorization Number (RMA) is obtained prior to shipment. Please call (1) 952-895-8292 or email: ctssales@ctshplc.com to get an RMA.